



## Which can be used on what?

	0 - 3 years			3 - 6 years					6 - 12 years				A							
	Panda size 1	Panda size 2	Wombat size 1	Serval size 1	Panther size 1	Panther size 2	Panda size 2½	Panda size 3	Wombat size 2	Serval size 2	Cougar 1	Cougar 2	Cougar 3	Panther size 2	Panther size 3	Panther size 4	Panda size 3	Panda size 4	Wombat size 3	Adult
Vests size 0	•			•																
Vests size ¾		•		•																
Vests size 1				•	•	•	•		•					•						
Vests size 2								•	•	•				•	•	•	•	•		
Vests size 3												•	•			•				
Vests size 4																				•
Cross vest size 0	•	•		•																
Cross vest size 1		•		•	•	•	•							•						
Cross vest size 2						•		•	•	•				•	•	•	•	•		
Cross vest size 3												•	•			•				
Cross vest size 4																				•
Belt Harness size 1	•	•	•	•	•	•	•													
Belt Harness size 2		•	•	•	•	•	•	•	•					•	•					
Belt Harness size 3							•	•	•	•				•	•	•	•			
Neoprene vest size 1	•	•		•																
Neoprene vest size 2		•		•	•	•	•		•	•				•						
Neoprene vest size 3						•		•	•	•	•	•	•	•	•	•	•	•		•
Sternum belt size 0	•	•		•																
Sternum belt size 1		•		•	•		•		•	•				•						
Sternum belt size 2						•		•	•	•	•	•	•	•	•	•	•	•		•
Chest belt, upholstered size 1	•	•	•		•	•	•							•						
Chest belt, upholstered size 2								•	•	•					•		•	•		
Chest belt, upholstered size 3											•			•	•		•	•		
Chest belt, upholstered size 4											•	•							•	•
5-point belt size 1					•	•	•		•											
5-point belt size 2							•	•						•			•			
Pelvis harness size 0	•	•		•																
Pelvis harness size 1	•	•		•																
Pelvis harness size 2					•	•	•	•		•				•			•			
Pelvis harness size 3															•	•		•		•
Y-hip strap, upholstered size 1	•	•		•																
Y-hip strap, upholstered size 2			•	•	•	•	•	•	•					•			•			
Y-hip strap, upholstered size 3									•	•					•	•		•	•	
Y-hip strap, upholstered size 4											•	•							•	•
Belt, 25 mm size 2	•	•	•	•																
Belt, 25 mm size 3			•	•	•	•	•	•	•					•			•			
Belt, 25 mm size 4									•	•					•	•		•		
Belt for transportation size 1					•															
Belt for transportation size 2						•								•						
Belt for transportation size 3															•	•				