

R82 Stingray - E1233

USD Retail Price List/ Order Form



Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

US Configuration


Effective: 02/01/2019 - Rev 7

- Frame (Stingray, Trekker, High-Low:v)
- Stingray Seat
- Cushions
- Head Support
- Side Support
- Hip Support
- Pommel
- Vests and Belts
- Tray
- Transport Fittings



- A (Back Height) _____
- B (Arm Height) _____
- C (Lower Leg Length) _____
- D (Seat Depth) _____
- E (Chest Width) _____
- F (Seat Width) _____
- G (Height: Seat to Top of Head) _____
- H (Shoulder Width) _____
- Age _____
- Weight _____

Measurements



	Size 1 Inches	Size 2 Inches
Effective Seat Width (A)	7-11 ³ / ₄ "	9-13 ³ / ₄ "
Seat Width Between Knee Supports	11"	12 ³ / ₄ "
Seat Width Between Hip Supports (Accessories)	7", 8 ¹ / ₂ ", 10 ¹ / ₄ "	9", 10 ¹ / ₂ ", 12"
Effective Seat Depth (B)	7 - 11 ³ / ₄ "	9 ³ / ₄ - 15 ¹ / ₂ "
Back Height at Shoulders (C)	10 ¹ / ₂ - 16"	14 - 19 ³ / ₄ "
Backrest Height (C1)	16 - 21 ¹ / ₂ "	19 ¹ / ₂ - 25 ¹ / ₄ "
Footrest to Seat Distance (D)	5 ¹ / ₂ - 17 ¹ / ₂ "	5 ¹ / ₂ - 17 ¹ / ₂ "
Seat Surface Height at Front Edge (E)	20 ¹ / ₄ "	20 ¹ / ₄ "
Back Width, Inside (F)	11 ³ / ₄ "	14"
Back Width, Outside (F1)	17 ¹ / ₂ "	20 ¹ / ₄ "
Overall Width (G)	23"	23"
Length with Turnable Front Wheels (H1)	28 ³ / ₄ "	28 ³ / ₄ "
Length with Large Front Wheels	32 ³ / ₄ "	32 ³ / ₄ "
Overall Height, Push Brace 15° (I)	39", 40", 41 ¹ / ₄ "	39", 40", 41 ¹ / ₄ "
Foot Plate, Width x Depth	12" x 7 ³ / ₄ "	12" x 7 ³ / ₄ "
Folded Length	35"	35"
Folded Height	19"	19"
Total Mass, Frame	29 ¹ / ₂ lbs.	29 ¹ / ₂ lbs.
Total Mass, Seat	12 lbs.	14 ¹ / ₂ lbs.
Max Load	88 lbs.	88 lbs.
Back Recline / Seat Plane Angle	Back: 170° / Seat: 45°	Back: 170° / Seat: 45°
Approx. Minimum User Height	29 ¹ / ₂ "	35 ¹ / ₂ "
Approx. Maximum User Height	50"	57 ¹ / ₂ "


NOTES:

Stingray Standard Features

- Tilt-in-Space - 45°
- Back Recline - 170°
- Turnable Seat - 180°
- One-Size Frame
- Accepts Two Seat Sizes for Maximum Growth Potential
- Adjustable Seat Depth
- Adjustable Back Height
- Suspension



Description	Part No.	USD Price	HPCPS
<input type="checkbox"/> Solid Tires, Individual Wheel Locks	8724212	\$2,445	E1233
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Hand Brakes ••	8724211	\$2,820	E1233
<input type="checkbox"/> Solid Tires, Central Wheel Locks	8724222	\$2,545	E1233
<input type="checkbox"/> Solid Tires, Central Wheel Locks & Hand Brakes	8724221	\$2,920	E1233
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Quick Release Rear Axels ••Custom order; Approx. lead time: 6-8 weeks. Non-returnable item	8722212	\$2,695	E1233
<input type="checkbox"/> Transit Option* *Required if using O2 Tank Holder	One Size 8720902	\$200	
<input type="checkbox"/> Medical Essential Carrier	One Size 8720600	\$97	
<input type="checkbox"/> Vent Tray ♦* 13" x 17" ♦Requires Central Wheel Locks *44 lbs. weight capacity	One Size 8776290	\$460	E1029
<input type="checkbox"/> O2 Tank Holder * •• *Transit Option (8720902) is required ••Custom order; Approx. lead time: 6-8 weeks. Non-returnable item	8776295	\$350	E2208
<input type="checkbox"/> IV Pole (height adj. & collapsible) (TR) Anti-Tipper Required	One Size C904566	\$200	K0105
<input type="checkbox"/> IV Pole Holder** **Required if using IV Pole	8776291	\$140	




Trekker Frame Standard Features


- One Piece Height Adjustable Push Handle (Inside Mount)
- Tilt-in-Space (-10° to +45°)
- 180° Reversible Seating Configuration
- Single Foot Activated Wheel Lock Bar
- 7½" x 2" Front and 11½" x 2½" Rear Solid Tires
- Quick Release Rear Wheels
- Made in the USA
- Transit Option - Transit Ready Frame Only



Description	Part No.	USD Price	HPCPS
<input type="checkbox"/> Trekker Frame - Transit Ready (TR)^^	C905362	\$2,595	E1234
<input type="checkbox"/> Stingray Seat Module - Interface Adaptor Kit for Trekker Base ♦ ♦Required with Trekker Frame ^ ^Stingray Seat Size 1 only compatible with Trekker Frame	C905363	N/C	
<input type="checkbox"/> Using other manufacturer's Head Support ♦Requires Head Support & Pelvic Positioning Belt	C995681	N/C	
<input type="checkbox"/> Height Adjustable Push Handle w/Outside Mount (TR) ^Trekker Frame comes standard with push buttons mounted inside Height Adjustable Push Handle	C904688	N/C	
<input type="checkbox"/> 12½ x 2" Rear Solid Knobby Tire (TR)••	C904117	\$199	
<input type="checkbox"/> 7½ x 2" Front & 11½ x 2½" Rear Pneumatic Knobby Tires (TR)•• ••Custom order; Approx. lead time: 6-8 weeks	C904118	\$299	

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> IV Pole (height adj. & collapsible) (TR) Anti-Tipper Required		One Size	C904566	\$200	K0105
	<input type="checkbox"/> O2 Tank Holder (TR) *°		One Size	C904565	\$175	E2208
	<input type="checkbox"/> LTV Ventilator Holder (TR) °		One Size	C904671	\$450	E1029
	<input type="checkbox"/> Ventilator/Suction/Hard Tray (TR)		One Size	C904595	\$450	E1029
	<input type="checkbox"/> Extended Ventilator/Suction/Hard Tray (TR)		One Size	C905326	\$525	E1029
	<input type="checkbox"/> Curb Tipper (TR)		One Size	C905310	\$98	
	<input type="checkbox"/> Rear Anti-Tipper (TR)		One Size	C904649	\$115	E0971
	<input type="checkbox"/> Under Seat Storage Basket (TR)		One Size	C904605	\$120	
	<input type="checkbox"/> Medical Necessity Storage Basket (TR)		One Size	C904606	\$295	
	* Must be removed during transit					
	° Manufacturer Installation Required					

High-Low Frame (Trekker)	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> High-Low:x sz 1, Manual Lift^		1,2	8911603-02	\$1,580	
	<input type="checkbox"/> Push Brace - Fixed Angle		One Size	8910545	\$160	
	<input type="checkbox"/> Push Brace - Angle Adjustable		1,2	8910555	\$385	
	^ Interface Adaptor Kit Required					

High-Low Frame	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> High-Low:v Indoor Frame		1	8911871-20	\$1,770	
	<input type="checkbox"/> High-Low:v Indoor Frame		1, 2	8911872-20	\$2,193	
	<input type="checkbox"/> Push Brace - Fixed Angle		One Size	8910545	\$160	
	<input type="checkbox"/> Push Brace, Angle Adjustable		One Size	8910555	\$385	

Stingray Seat	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Seat and Back size 1^^		1	8711120		
	<input type="checkbox"/> Seat size 1		1	x8711000	\$410	E2292
	<input type="checkbox"/> Back size 1		1	x8712000	\$505	E2291
	<input type="checkbox"/> Seat and Back size 2		2	8712220		
	<input type="checkbox"/> Seat size 2		2	x8721000	\$410	E2292
	<input type="checkbox"/> Back size 2		2	x8722000	\$505	E2291
	<input checked="" type="checkbox"/> Fully Reclining Back (Required)		One Size	x8702000	\$405	E1014

Adjustable seat depth and height adjustable back
 ^^ **Stingray Seat Size 1 only compatible with Trekker Frame**

Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------

Required, must select a color

X = Color

- Blue = 1 Red = 2
 Pink = 5 Grey = 4

- Seat and Back Cushion (Pr.) 1 875051x N/C
 Seat and Back Cushion (Pr.) 2 875052x N/C










* **Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item**

Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------

- Footrest Rod, Elevating ♦ 8¼" Short x8711100 \$110 K0108
 Footrest Rod, Elevating ♦ 11" Long x8712200 \$110 K0108
 ♦ **Required with Angle Adjustable Footplate (x8710010)**

- Angle Adjustable Footplate ♦ One Size x8710010 \$95 K0040
 ♦ **Requires Footrest Rod Short or Long (x8711100 or x8712200)**



Accessories	Description	Measurements	Size	Part No.	USD Price	HPCPS
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.) ♦		1	8740250	\$225	E0956
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) (Required) ♦		1	x8740250	\$225	E1028
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.) ♦		2	8740260	\$225	E0956
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) (Required) ♦		2	x8740260	\$225	E1028
♦ Hip Support and Multi-Adjustable Hardware must be ordered together						
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	4" x 4"	1	8720301	\$195	E0956
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	5" x 5¾"	2	8720302	\$195	E0956
<hr/>						
Description		Measurements	Size	Part No.	USD Price	HPCPS
<input type="checkbox"/> Covers for Knobs on Side Supports (Pr.) ●●			One Size	8720360	\$30	
●● Custom order; Approx. lead time: 6-8 weeks. Non-returnable item						
	<input type="checkbox"/> Medial Thigh Support (Ea.) ♦	3½" x 1½" x 2¾"	1	8720501	\$135	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		1	x8720501	\$165	E1028
	<input type="checkbox"/> Medial Thigh Support (Ea.) **	3½" x 3½" x 3½"	2	8720502	\$135	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		2	x8720502	\$165	E1028
	<input type="checkbox"/> Medial Thigh Support (Ea.) **	3½" x 3½" x 4¾"	3	8720503	\$135	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		3	x8720503	\$165	E1028
♦ Medial Thigh Support and Adjustable, Removable Hardware must be ordered together						
	<input type="checkbox"/> Anti-Thrust Wedge, Black (Pr.)	4" x 6"	1, 2, 3	8910400-3	\$46	K0108
<hr/>						
Description		Measurements	Size	Part No.	USD Price	HPCPS
	<input type="checkbox"/> Contoured Headrest, Black		1	8740810	\$200	E0955
	<input type="checkbox"/> Contoured Headrest, Black		2	8740820	\$200	E0955
	<input type="checkbox"/> Sculpted Head Support ①	9½" x 4½"	Small	8720831	\$205	E0955
	<input type="checkbox"/> Sculpted Head Support ①	11" x 5¾"	Large	8720832	\$205	E0955
	① Preferred option for Stingray seating					
	<input type="checkbox"/> Anatomic Head Support ♦	A:5½" x B:5½" x C:9½"	Small	5099001	\$255	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C	
	<input type="checkbox"/> Anatomic Head Support ♦	A:6" x B:5½" x C:10"	Medium	5099002	\$255	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C	
	<input type="checkbox"/> Anatomic Head Support ♦	A:6½" x B:5½" x C:12½"	Large	5099003	\$255	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C	
♦ Turtle Bar System Short or Long Required; 98619-25 or 98619-26						
Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable						
	<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$230	E1028
	<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$230	E1028
♦ Required with Anatomic Head Support; 5099001, 5099002 or 5099003						
	<input type="checkbox"/> Headrest Fittings ♦		One Size	8720850	\$51	K0108
♦ Required with Turtle Bar System or Universal Headrest Mount						

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	-------------	--------------	------	----------	-----------	-------



2-Part Headrest
(Mounts directly onto backrest) One Size 8720870 \$270 E0955

Universal Headrest Mount One Size 87590 \$61 E1028


◆ Required when using other manufactures Head Support

Allow most Whitmyer, Stealth, Otto Bock headrest, I2I to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation


Other types of headrest available, please see Head Supports Order Form

Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	--------------	------	----------	-----------	-------


■ XX = 60 (Black)



Cross Vest (Ea.) A: 4 1/4 x B: 1/2 x C: 11" 0 89625-0XX \$190 E0960
 Cross Vest (Ea.) A: 5 x B: 1/2 x C: 12" 1 89625-1XX \$190 E0960
 Cross Vest (Ea.) A: 5 1/2 x B: 1/2 x C: 13" 2 89625-2XX \$190 E0960
 Cross Vest (Ea.) A: 6 1/4 x B: 2 1/2 x C: 15 1/4" 3 89625-3XX \$190 E0960
 Cross Vest (Ea.) A: 7 1/2 x B: 2 1/2 x C: 18 1/2" 4 89625-4XX \$190 E0960




Vest (Ea.) 1 1/4" x 2 1/2" x 5 1/2" x 3 1/2" x 11" 0 89111-0XX \$200 E0960
 Vest (Ea.) 1 1/4" x 2 1/2" x 5 3/4" x 4 1/4" x 12" 3/4 89111-075XX \$200 E0960
 Vest (Ea.) 2 1/4" x 2 1/2" x 7" x 5" x 13 1/4" 1 89111-1XX \$200 E0960
 Vest (Ea.) 2 1/2" x 3 1/2" x 9 1/4" x 5 1/2" x 15 1/4" 2 89111-2XX \$200 E0960
 Vest (Ea.) 2 3/4" x 4 3/4" x 11 x 5 3/4" x 18 1/4" 3 89111-3XX \$200 E0960
 Vest (Ea.) 2 3/4" x 6 3/4" x 13 3/4" x 7 1/2" x 19" 4 89111-4XX \$200 E0960



Neoprene Vest (Ea.) A: 7 x B: 5 1/2 x C: 4 x D: 4 1/2" 1 89113-1 \$200 E0960
 Neoprene Vest (Ea.) A: 10 1/2 x B: 11 x C: 6 1/2 x D: 5 3/4" 2 89113-2 \$200 E0960
 Neoprene Vest (Ea.) A: 11 1/4 x B: 9 x C: 7 x D: 6 1/2" 3 89113-3 \$200 E0960




4-Point Padded Hip Belt L: 10 1/4-11 1/2" 1 85427-1 \$135 E0978
 4-Point Padded Hip Belt L: 12 1/2-13 1/2" 2 85427-2 \$135 E0978



2-Point Pelvic Belt, Black L: 13-14 1/2" 1 85424-1 \$97 E0978



Fix-Locks, Black, Bag of 6 Pcs ◆ One Size 862109-35 \$42
 ◆ Required for use of Vests



Foot Straps, Plastic Lock (Pr.) Max: 22", A: 13 1/2-20 1/2" One Size 85096 \$36 K0040
 Foot Straps, Velcro (Pr.) One Size 85102 \$51 K0040

Other types of vests and belts available, please see Vest and Belts Order Form

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Folding Top (Canopy), Black		1	8740110	\$330	
	<input type="checkbox"/> Folding Top (Canopy), Black		2	8740120	\$330	
	<input type="checkbox"/> Foam Back Plate ●●	5" x 12"	1	8776325	\$23	
	<input type="checkbox"/> Foam Back Plate ●●	5½" x 13¾"	2	8776326	\$23	
	Future Growth - Fills out the gap between seat and the back					
	●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item					
	<input type="checkbox"/> Tray with Edge, Grey ♦	17.3" x 13.3"	One Size	8776350	\$100	E0950
	<input type="checkbox"/> Adjustable Removable Tray Hardware ♦		One Size	x8776350	\$165	E1028
	♦ Tray and Adjustable Removable Hardware must be ordered together					
	<input type="checkbox"/> Rain Cover ♦		1	8720141	\$56	
	<input type="checkbox"/> Rain Cover ♦		2	8720142	\$56	
	♦ Requires Folding Top (Canopy)					
	<input type="checkbox"/> Handle Grip	14"	1	8720401	\$200	
	<input type="checkbox"/> Handle Grip	16"	2	8720402	\$200	
	<input type="checkbox"/> Backpack, Gray ●●		1	8720711	\$56	
	<input type="checkbox"/> Backpack, Gray ●●		2	8720721	\$56	
	●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item					

Spare Parts	Description	Measurements	Size	Part No.	USD Price	HCPCS
	Required, must select a color					
	X = Color					
	<input type="checkbox"/> Blue = 1	<input type="checkbox"/> Red = 2				
	<input type="checkbox"/> Pink = 5	<input type="checkbox"/> Grey = 4				
	<input type="checkbox"/> Seat and Back Cushion (Pr.) ●●		1	875051x	\$433	
	<input type="checkbox"/> Seat and Back Cushion (Pr.) ●●		2	875052x	\$433	
	●●Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item					
	<input type="checkbox"/> Adductor Covers, Black (Pr.)		1	8750410	\$56	
	<input type="checkbox"/> Adductor Covers, Black (Pr.)		2	8750420	\$56	
	<input type="checkbox"/> Hardware Kit for Hip Supports, Fixed		One Size	872023-KIT	\$35	

Spare Parts	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Release Strap for Back Angling		1	87102250	\$46	
	<input type="checkbox"/> Release Strap for Back Angling		2	87102255	\$46	
	<input type="checkbox"/> Footrest Rod, Elevating	8¼"	Short	8775317	\$110	
	<input type="checkbox"/> Footrest Rod, Elevating	11"	Long	8775318	\$110	
	<input type="checkbox"/> Angle Adjustable Footplate		One Size	8775010	\$260	

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.

The R82 Stingray mobility base was tested with a footrest and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.

For shipping rates, please call 1-844-876-6245
12801 E Independence Blvd, Matthews, NC 28105 • 1-844-876-6245 • Web: www.r82.com
Prices in USD & is subject to change without notice
For Quotes and Orders please send completed Order Form to: us.info@r82.com