

R82 Mustang - E8002

USD Retail Price List/ Order Form

Account #: _____ Purchase Order No.: _____ Date: _____


Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

US Configuration



Effective: 7.20.20 - Rev 6

- Mustang
- Accessories**
- Chest Support
- Head Support
- Back Support
- Arm Support
- Seat/Sling Seat
- Leg Separator
- Ankle Prompts



- A (Total Height) _____
- B (Shoulder Height) _____
- C (Chest Height) _____
- D (Buttock Height) _____
- E (Chest Width) _____
- F (Shoulder Width) _____
- Age _____
- Weight _____


Measurements



	Size 1 Inches	Size 2 Inches	Size 3 Inches	Size 4 Inches
Width, total (G)	20¼"	22½"	25½"	27½"
Width, inside frame (G1)	14¾"	17½"	19½"	22"
Length with 4" wheels (H)	29¼"	32½"	-	-
Length with 8" wheels (H)	32"	35"	39"	43"
Height with 4" wheels (I)	14¾-25¼"	23½-33½"	-	-
Height with 8" wheels (I)	19-25¼"	23½-33½"	27½-41"	34-50¾"
Seat*) height above floor, 4" wheels	6½-13½"	9½-16"	-	-
Seat*) height above floor, 8" wheels	10-17¼"	12¾-19½"	17½-26½"	24-32¾"
Width between handles *)	5½-13½"	5½-13½"	5½-23½"	15½-23½"
Angle range of center bar (a1)	30°	30°	0-25°	0-25°
Turning radius	18	21½	23½	25¼"
User height, floor to axil	14¾-25¼"	23½-33½"	27½-41½"	34-50¾"
Weight **) with 4" wheels	22 lbs.	24 lbs.	-	-
Weight **) with 8" wheels	28½ lbs.	31 lbs.	40 lbs.	42 lbs.
Max. load / user weight	66 lbs.	88 lbs.	132 lbs.	176 lbs.
Approx. min. user height (4" wheels)	26½"	33"	-	-
Approx. min. user height (8" wheels)	28¾"	34½"	40½"	49½"
Approx. max. user height (4" wheels)	36½"	44½"	-	-
Approx. max. user height (8" wheels)	36½"	44½"	54½"	66"


The Mustang is measured as standard delivered and with the center bar in vertical position. Different settings may result in another set of measurements. To accomplish some of the measurements above, a slight adjustment might be required. We supply various types of accessories, which is why the weight is listed without.

Mustang Standard Features



- Four Sizes with Four Solid Wheels
- Drag Wheel Function
- Reverse Brake On/Off
- Directional Stabilizer
- Angle Adjustment
- Height Adjustment
- Chest Support, Height/Width Adjustable

Description	Size	Part No.	USD Price	HPCPS
<input type="checkbox"/> Mustang, Yellow, 4" Casters	1	869011	\$1,935	E8002
<input type="checkbox"/> Mustang, Yellow, 8" Casters	1	869015	\$1,935	E8002
<input type="checkbox"/> Mustang, Yellow, 8" Casters Slim Line Wheels	1	869015-SL	\$1,935	E8002
<input type="checkbox"/> Mustang, White, 4" Casters	2	869025-4	\$2,015	E8002
<input type="checkbox"/> Mustang, White, 8" Casters	2	869024	\$2,015	E8002
<input type="checkbox"/> Mustang, White, 8" Casters Slim Line Wheels	2	869024-SL	\$2,015	E8002
<input type="checkbox"/> Mustang, Red, 8" Casters	3	869033	\$2,850	E8002
<input type="checkbox"/> Mustang, Red, 8" Casters Slim Line Wheels	3	869033-SL	\$2,850	E8002
<input type="checkbox"/> Mustang, Anthracite (Grey), 8" Casters	4	869042	\$3,260	E8002
<input type="checkbox"/> Mustang, Anthracite (Grey), 8" Casters Slim Line Wheels	4	869042-SL	\$3,260	E8002



Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------

Required, must select one:



	Inner Circum.	Width Range	Size	Part No.	USD Price
<input type="checkbox"/> Chest Support	15½-25¼"	6 - 8½"	Small	869810	N/C
<input type="checkbox"/> Chest Support	26¼ -30¾"	9½ - 12"	Medium	869820	N/C
<input type="checkbox"/> Chest Support	34¼ -44¾"	11 - 14¼"	Large	869830	N/C


Description	Size	Part No.	USD Price
-------------	------	----------	-----------



<input type="checkbox"/> Cover for Chest Support	H:3¼ x W:9¾"	Small	869810-160	\$115
<input type="checkbox"/> Cover for Chest Support	H:4¾ x W:14¾"	Medium	869820-260	\$140
<input type="checkbox"/> Cover for Chest Support ●●	H:5¾ x W:19"	Large	869830-360	\$150

●●Custom order; Approx. lead time: 6-8 weeks (sz Large)

Accessories	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



	Inner Circum.	Width Range	Size	Part No.	USD Price
<input type="checkbox"/> Pelvic Support	15½-25¼"	6 - 8½"	Small	869610	\$475
<input type="checkbox"/> Pelvic Support	26¼ -30¾"	9½ - 12"	Medium	869620	\$490
<input type="checkbox"/> Pelvic Support ●●	34¼ -44¾"	11 - 14¼"	Large	869630	\$505

●●Custom order; Approx. lead time: 6-8 weeks (sz Large)




<input type="checkbox"/> Extra Pad for Chest/Pelvic Sup.	H:4¾ x W:2¼"	One Size	869225	\$61
--	--------------	----------	--------	------





Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------

<input type="checkbox"/> Sling Seat for Small Chest Support	W:7 x L:9½"	1	869570	\$175
<input type="checkbox"/> Sling Seat for Med Chest Support	W:7 x L:12¾"	2	869571	\$175






Note: Sling Seat n/a on size's 3 and 4 Mustang



<input type="checkbox"/> Mustang Seat	H:6¾ x D:10"	1, 2	869110	\$255
<input type="checkbox"/> Mustang Seat	H:7"x D:12¾"	3, 4	869120	\$270
<input type="checkbox"/> Mustang Comfort Seat	L:8½"	1, 2	869115	\$463
<input type="checkbox"/> Mustang Comfort Seat	L:11¾"	3, 4	869125	\$480



Accessories	Description	Measurements	Size	Part No.	USD Price
	<input type="checkbox"/> Hip Support on Curved Bar ♦		1, 2	869222	\$385
	<input type="checkbox"/> Hip Support on Curved Bar ♦		3, 4	869232	\$385
	♦Requires Mustang Seat or Mustang Comfort Seat				
	<input type="checkbox"/> Hip Sup. w/ Bracket for Head Sup ♦ W:7¼"-10"		1, 2	869221	\$385
	<input type="checkbox"/> Hip Sup. w/ Bracket for Head Sup ♦ W:10-12¾"		3, 4	869231	\$385
	♦Requires Mustang Seat or Mustang Comfort Seat				
	<input type="checkbox"/> Sacrum Support ♦♦	H:7", D:5-8¼"	Small	869250	\$160
	<input type="checkbox"/> Sacrum Support ♦♦	H:8", D:6½"-9½"	Large	869255	\$220
	♦Requires Mustang Seat or Mustang Comfort Seat				
	♦♦Custom order; Approx. lead time: 6-8 weeks (sz Large)				
	Note: Measurements from Seat				
	<input type="checkbox"/> Abdominal Support	H:5 x W:7½"	One Size	869260	\$92





Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------


	<input type="checkbox"/> Hand Bar for Anterior Use	H:41/4 x W:8¼"	Small	869360	\$210
	<input type="checkbox"/> Hand Bar for Anterior Use ♦♦	H:41/4 x W:5½"	Large	869365	\$210
	♦♦Custom order; Approx. lead time: 6-8 weeks (sz Large)				
	<input type="checkbox"/> Armrest for Anterior Use w/ Handles		Small	869350	\$410
	<input type="checkbox"/> Armrest for Anterior Use w/ Handles		Large	869355	\$505
	<input type="checkbox"/> Arm Prompts for Anterior Use		Small	869310	\$385
	<input type="checkbox"/> Arm Prompts for Anterior Use		Medium	869320	\$385
	<input type="checkbox"/> Arm Prompts for Anterior Use		Large	869325	\$405
	<input type="checkbox"/> Arm Prompts for Posterior Use		Small	869330	\$405
	<input type="checkbox"/> Arm Prompts for Posterior Use ♦♦		Medium	869340	\$405
	<input type="checkbox"/> Arm Prompts for Posterior Use ♦♦		Large	869345	\$405
	♦♦Custom order; Approx. lead time: 6-8 weeks (sz Medium and Large)				
	<input type="checkbox"/> Padded Forearm Positioner (Ea.)		One Size	869370	\$68

QTY

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------

	<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$230
	<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$230
	♦Required with Anatomic Head Support				
	<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$255
	<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C
	<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$255
	<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C
	<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$255
	<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C
	♦Turtle Bar System Short or Long Required; 98619-25 or 98619-26				
	<input type="checkbox"/> Headband System ♦		One Size	99500	\$230
	♦For Anatomic Head Support Only				

Accessories	Description	Measurements	Size	Part No.	USD Price																				
	<input type="checkbox"/> Connecting Rod Between Hip and Head Support ♦ ♦Required to connect Headrest to Hip or Sacrum Supports		One Size	869226	\$56																				
	<input type="checkbox"/> Leg Separator for sz 1 w/ 4" Wheels		1	869410	\$290																				
	<input type="checkbox"/> Leg Separator for sz 1 w/ 8" Wheels		1	869410-8	\$290																				
	<input type="checkbox"/> Leg Separator for sz 2 w/ 4" Wheels		2	869430	\$290																				
	<input type="checkbox"/> Leg Separator for sz 2 w/ 8" Wheels		2	869420	\$290																				
	<input type="checkbox"/> Leg Separator for sz 3		3	869430	\$290																				
	<input type="checkbox"/> Leg Separator for sz 4		4	869440	\$350																				
	<input type="checkbox"/> Ankle Prompts for sz 1 ♦		1	869450	\$175																				
	<input type="checkbox"/> Ankle Prompts for sz 2 w/ 8" Wheels		2	869455	\$175																				
	<input type="checkbox"/> Ankle Prompts for sz 2 w/ 4" Wheels		2	869455-4	\$175																				
	<input type="checkbox"/> Ankle Prompts for sz 3		3	869460	\$190																				
	<input type="checkbox"/> Ankle Prompts for sz 4		4	869465	\$190																				
		♦Can be used on Mustang sz 1 with both 4" and 8" Wheels																							
	<input type="checkbox"/> Leg Pad		One Size	869470	\$135																				
<table border="1"> <thead> <tr> <th>Description</th> <th>Measurements</th> <th>Size</th> <th>Part No.</th> <th>USD Price</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/> Attendant Guide Bar <i>Note: N/A on Mustang sz 3 and 4</i></td> <td>L: 27¼"</td> <td>1, 2</td> <td>869510</td> <td>\$200</td> </tr> <tr> <td><input type="checkbox"/> Storage Pouch ●● ●●Custom order; Approx. lead time: 6-8 weeks</td> <td>H: 5¾ x W: 5"</td> <td>1, 2</td> <td>869560</td> <td>\$97</td> </tr> <tr> <td><input type="checkbox"/> Storage Bag ●● ●●Custom order; Approx. lead time: 6-8 weeks</td> <td>H: 12 x W: 15½"</td> <td>3, 4</td> <td>869565</td> <td>\$195</td> </tr> </tbody> </table>						Description	Measurements	Size	Part No.	USD Price	<input type="checkbox"/> Attendant Guide Bar <i>Note: N/A on Mustang sz 3 and 4</i>	L: 27¼"	1, 2	869510	\$200	<input type="checkbox"/> Storage Pouch ●● ●●Custom order; Approx. lead time: 6-8 weeks	H: 5¾ x W: 5"	1, 2	869560	\$97	<input type="checkbox"/> Storage Bag ●● ●●Custom order; Approx. lead time: 6-8 weeks	H: 12 x W: 15½"	3, 4	869565	\$195
Description	Measurements	Size	Part No.	USD Price																					
<input type="checkbox"/> Attendant Guide Bar <i>Note: N/A on Mustang sz 3 and 4</i>	L: 27¼"	1, 2	869510	\$200																					
<input type="checkbox"/> Storage Pouch ●● ●●Custom order; Approx. lead time: 6-8 weeks	H: 5¾ x W: 5"	1, 2	869560	\$97																					
<input type="checkbox"/> Storage Bag ●● ●●Custom order; Approx. lead time: 6-8 weeks	H: 12 x W: 15½"	3, 4	869565	\$195																					

Spare Parts	Description	Measurements	Size	Part No.	USD Price
	<input type="checkbox"/> 4" Wheel and Holder, Left Side		1	79436-L-KIT	\$285
	<input type="checkbox"/> 4" Wheel and Holder, Right Side		1	79436-R-KIT	\$285
	<input type="checkbox"/> 8" Wheel and Holder, Left Side		2, 3, 4	79438-L-KIT	\$410
	<input type="checkbox"/> 8" Wheel and Holder, Right Side		2, 3, 4	79438-R-KIT	\$410
	<input type="checkbox"/> 8" Slim-Line Wheel and Holder, Left Side		2, 3, 4	79439-L-KIT	\$410
	<input type="checkbox"/> 8" Slim-Line Wheel and Holder, Right Side		2, 3, 4	79439-R-KIT	\$410
	<input type="checkbox"/> 8" Slim-Line Wheels and Bushings, Kit of 4		2, 3, 4	79439-SLKIT	\$408

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in USD & is subject to change without notice

For Quotes and Orders please send completed Order Form to: sales.us@etac.com or fax: 1-800-760-0537