

# R82 Wombat Living

## USD Retail Price List/ Order Form

Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

**US Configuration**


**Effective: 04.27.20 - Rev 7**

- Wombat Living
- Head Support
- Side Support
- Armrest
- Hip Support
- Pommel
- Vests and Belts
- Tray
- Helper's Handle




- A** (Back Height) \_\_\_\_\_
- B** (Arm Height) \_\_\_\_\_
- C** (Lower Leg Length) \_\_\_\_\_
- D** (Seat Depth) \_\_\_\_\_
- E** (Chest Width) \_\_\_\_\_
- F** (Seat Width) \_\_\_\_\_
- G** (Height from Seat to Top of Head) \_\_\_\_\_
- H** (Shoulder Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

**Measurements**



**Size 1**  
Inches



**Size 2**  
Inches

**Size 3**  
Inches

Effective seat width <b>(A)</b>	11¾"	14½"	16½"
Seat width between hip supports	4¼-11¾"	6¾-14½"	8½-16½"
Seat width, pelvis support w/sides	6¾-9½"	9-11¾"	11½-14¼"
Width, pelvis support without sides	9½"	9½"	11½"
Effective seat depth <b>(B)</b>	7-12¼"	9¾-15"	12½-19¾"
Backrest height <b>(C)</b>	11-13½"	12¼-15"	15¼-19¾"
Backrest height <b>(C)</b> (high back)	14¼-17"	15½-18"	18-22¾"
Footrest to seat distance <b>(D)</b>	6-10"	6-13¾"	6-17¼"
Seat surface height at front edge (4" Castors)	15¼-26¾"	15¼-26¾"	16½-27½"
Back width, flat cushion <b>(F)</b>	9¾"	11¾"	13¾"
Back width, contoured cushion, between side supports <b>(F)</b>	6¾"	9"	10½"
Arm rest height above seat cushion long rod and pad	2½-7"	2½-7"	2½-7"
Backrest angle <b>(a1)</b>	2½-10½"	2½-10½"	2½-10½"
Seat plane angle <b>(a2)</b>	90° to 130°	90° to 130°	90° to 133°
Leg to seat surface angle <b>(a3)</b>	-15° to 20°	-15° to 20°	-10° to 22°
Overall width <b>(G)</b>	-30° to 30°	-30° to 30°	-30° to 30°
Length <b>(H)</b>	21½"	21½"	24"
Max load	23½"	23½"	28¼"
Total mass	99.2 lbs.	99.2 lbs.	154.3 lbs.
	34.1 lbs.	34.1 lbs.	48½ lbs.

**NOTES:**

**Wombat Living Standard Features**



- Tilt in Space
- Back Recline and Prone Angle
- Adjustable Seat Depth
- Adjustable Back Height
- Height and Angle Adjustable Foot Support
- Without Seat/Back Cushions

**Description Pricing**

			95XXXX-00
<input type="checkbox"/> Wombat Living, size 1	\$3,125	1	
<input type="checkbox"/> Wombat Living, size 2	\$3,125	2	
<input type="checkbox"/> Wombat Living, size 3	\$3,125	3	
<input type="checkbox"/> Wombat Living, size 1 with High Back●●	\$3,125	5	
<input type="checkbox"/> Wombat Living, size 2 with High Back●●	\$3,125	6	
<input type="checkbox"/> Wombat Living, size 3 with High Back●●	\$3,125	7	
●● <i>Custom Order; Approx. lead time: 6-8 weeks</i>			
<input type="checkbox"/> Frame Color: Anthracite (Dark Grey)	\$0	2	
<input type="checkbox"/> 4" Front/Rear Wheels	\$0	2	
<input type="checkbox"/> 4" Front/Rear Wheels with central brake - <i>upcharge</i> ●●	\$235	3	
●● <i>Custom Order; Approx. lead time: 6-8 weeks</i>			
<input type="checkbox"/> w/Gas Spring: Height Pedal on Frame, Tilt Lever Under Seat	\$0	1	
<input type="checkbox"/> w/Gas Spring: Height Pedal on Frame, Tilt Lever on Push Brace - <i>upcharge</i>	\$265	2	
<input type="checkbox"/> w/power actuator: Powered height setting, tilt lever under seat - <i>upcharge</i> ●●	\$1,420	5	
<input type="checkbox"/> w/power actuator: Powered height setting, tilt lever on push brace - <i>upcharge</i> ●●	\$1,420	6	
●● <i>Custom Order for frame size 1 &amp; 2; Approx. lead time: 6-8 weeks</i>			

**Important!**

Not recommended for users with a lot of involuntary movements

**Accessories Description Size Part No. USD Price**

**Required, must select a color**

**xx = Color**


- Red = 24       Petroleum = 25       Grey = 28



<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	1	955031-2xx	\$140
<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	2	955032-2xx	\$140
<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	3	955033-2xx	\$140



<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	1	955031-238	\$240
<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	2	955032-238	\$240
<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	3	955033-238	\$240



**Description Size Part No. USD Price**

<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	1	955021-2xx	\$270
<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	2	955022-2xx	\$270
<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	3	955023-2xx	\$270


**Accessories Description Size Part No. USD Price**



<input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.)	1	955021-238	\$270
<input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.)	2	955022-238	\$270



Seat Cushion, Contoured w/Grey Cover (Ea.) 3 955023-238 \$270



Seat Cushion, Contoured w/ Incontinence Cover (Ea.) 1 955021-237 \$290  
 Seat Cushion, Contoured w/ Incontinence Cover (Ea.) 2 955022-237 \$290  
 Seat Cushion, Contoured w/ Incontinence Cover (Ea.) 3 955023-237 \$290

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



Lower Back Ext, Black (Ea.) \* H:2½" x W:10½" 3 95140-60 \$97  
*\*Only available on Wombat Living sz 3*



Lumbar Foam Wedge (Ea.) 1 955030-1 \$25  
 Lumbar Foam Wedge (Ea.) 2 955030-2 \$25  
 Lumbar Foam Wedge (Ea.) 3 955030-3 \$25

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------




Hip Support (Pr.) H:3¼" x W:7½" 1,2 95110-1 \$200  
 Covers for Hip Supports (Required) 1,2 95110-60 \$0  
 Hip Support (Pr.) H:4" x W:10" 3 95111-1 \$200  
 Covers for Hip Supports (Required) 3 95111-60 \$0



Pelvis Support without Sides (Ea.) ♦ 1,2 95140-61 \$231  
 Pelvis Support without Sides (Ea.) ♦ 3 95140-63 \$231  
*♦ Only available on Wombat Living with High Back*



Pelvis Support with Sides (Ea.) ♦ 1 95140-64 \$417  
 Pelvis Support with Sides (Ea.) ♦ 2 95140-65 \$417  
 Pelvis Support with Sides (Ea.) ♦ 3 95140-66 \$417  
*♦ Only available on Wombat Living with High Back*



Pommel (Ea.) H:3½" x W:3½" x D:3½" 1 95130-1 \$145  
 Pommel (Ea.) H:3½" x W:3½" x D:4½" 2 95130-2 \$145  
 Pommel (Ea.) H:3½" x W:3½" x D:6" 3 95130-3 \$145  
 Pommel (Ea.) H:3½" x W:3½" x D:7" 4 95130-4 \$145




Knee Abduction Pommel One Size 95230-10 \$170



Padding for Knee Abduction Pommel One Size 95430-10260 \$65

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



Armrest, Short (Pr.) 1,2 95112-1 \$350  
 Armrest, Short (Pr.) 3 95114 \$350






Armrest, Long Rod and Pad (Pr.) 1,2 95113-1 \$350  
 Armrest, Long Rod and Pad (Pr.) 3 95115 \$350

Accessories	Description	Size	Part No.	USD Price
-------------	-------------	------	----------	-----------




Fixed Side Supports w/Covers (Pr.) 1,2 95121 \$170

Swing-Away Lateral Support Bracket, RH (Ea.) One Size 83400-R \$225


<input type="checkbox"/> Swing-Away Lateral Support Bracket, LH (Ea.)	One Size	83400-L	\$225	
<input type="checkbox"/> Swing-Away Lateral Support Bracket, Offset, RH (Ea.)	One Size	83600-R	\$225	
<input type="checkbox"/> Swing-Away Lateral Support Bracket, Offset, LH (Ea.)	One Size	83600-L	\$225	
 <input type="checkbox"/> Bar, Straight (Ea.)	4 <sup>3</sup> / <sub>4</sub> "	One Size	81932	\$0
<input type="checkbox"/> Bar, Straight, Stainless (Ea.)	3 <sup>3</sup> / <sub>4</sub> "	One Size	82569	\$0
 <input type="checkbox"/> Bar, Offset, 1" (Ea.)	4 <sup>3</sup> / <sub>4</sub> "	One Size	81968	\$0
<input type="checkbox"/> Bar, Offset, 1 <sup>3</sup> / <sub>4</sub> " (Ea.)	4 <sup>3</sup> / <sub>4</sub> "	One Size	82000	\$78
<input type="checkbox"/> Bar, Offset, 2 <sup>1</sup> / <sub>2</sub> " (Ea.)	4 <sup>3</sup> / <sub>4</sub> "	One Size	81969	\$0
 <input type="checkbox"/> Curved Cushion (Ea.)	3 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>2</sub> "	One Size	6175406	\$105
<input type="checkbox"/> Curved Cushion (Ea.)	2 <sup>1</sup> / <sub>2</sub> " x 6 <sup>3</sup> / <sub>4</sub> "	One Size	6175405	\$105
<input type="checkbox"/> Curved Cushion (Ea.)	2 <sup>1</sup> / <sub>2</sub> " x 6"	One Size	6175404	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	6" x 4 <sup>3</sup> / <sub>4</sub> "	One Size	6175403	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	4" x 4"	One Size	6175402	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	2 <sup>2</sup> / <sub>3</sub> " x 4 <sup>1</sup> / <sub>2</sub> "	One Size	6175401	\$105

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------


■ XX = 60 (Black)




<input type="checkbox"/> Neoprene Vest (Ea.)	A:7 x B:5 <sup>1</sup> / <sub>2</sub> x C:4 x D:4 <sup>1</sup> / <sub>2</sub> "	1	89113-1	\$200
<input type="checkbox"/> Neoprene Vest (Ea.)	A:10 <sup>1</sup> / <sub>2</sub> x B:11 x C:6 <sup>1</sup> / <sub>2</sub> x D:5 <sup>3</sup> / <sub>4</sub> "	2	89113-2	\$200
<input type="checkbox"/> Neoprene Vest (Ea.)	A:11 <sup>1</sup> / <sub>4</sub> x B:9 x C:7 x D:6 <sup>1</sup> / <sub>2</sub> "	3	89113-3	\$200




<input type="checkbox"/> Cross Vest (Ea.)	A: 5 x B:1 <sup>1</sup> / <sub>2</sub> x C:12"	1	89625-1XX	\$190
<input type="checkbox"/> Cross Vest (Ea.)	A: 5 <sup>1</sup> / <sub>2</sub> x B:1 <sup>1</sup> / <sub>2</sub> x C:13"	2	89625-2XX	\$190
<input type="checkbox"/> Cross Vest (Ea.)	A: 6 <sup>1</sup> / <sub>4</sub> x B:2 <sup>1</sup> / <sub>2</sub> x C:15 <sup>1</sup> / <sub>4</sub> "	3	89625-3XX	\$190



<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:11 <sup>1</sup> / <sub>4</sub> x B:7 <sup>1</sup> / <sub>2</sub> x C:9 <sup>1</sup> / <sub>4</sub> -12"	1	85428-1	\$270
<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:14xB:10xC:10-12 <sup>3</sup> / <sub>4</sub> "	2	85428-2	\$270
<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:16 <sup>1</sup> / <sub>2</sub> x B:12 <sup>1</sup> / <sub>2</sub> x C:10-13"	3	85428-3	\$270



<input type="checkbox"/> 4-Point Padded Hip Belt	L:10 <sup>1</sup> / <sub>4</sub> -11 <sup>1</sup> / <sub>2</sub> "	1	85427-1	\$135
<input type="checkbox"/> 4-Point Padded Hip Belt	L: 12 <sup>1</sup> / <sub>2</sub> -13 <sup>1</sup> / <sub>2</sub> "	2	85427-2	\$135



<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14 <sup>1</sup> / <sub>2</sub> "	1	85424-1	\$97
---	--	---	---------	------




<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ** **Required for use of Vests	One Size	862109-35	\$42
---	----------	-----------	------

<input type="checkbox"/> Y-Hip Strap, Upholstered, 1" End Girth (Ea.)	1	85427-1	\$135
---	---	---------	-------

<input type="checkbox"/> Y-Hip Strap, Upholstered, 1" End Girth (Ea.)	2	85427-2	\$135
---	---	---------	-------

Accessories	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 5 <sup>1</sup> / <sub>2</sub> -7 x B:8 <sup>1</sup> / <sub>2</sub> -10"	1	85426-1	\$185
<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 6 <sup>1</sup> / <sub>4</sub> -7 <sup>3</sup> / <sub>4</sub> x B:10-11 <sup>3</sup> / <sub>4</sub> "	2	85426-2	\$185
<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 7-8 <sup>1</sup> / <sub>2</sub> x B:11-12 <sup>1</sup> / <sub>2</sub> "	3	85426-3	\$185





Neoprene Sandals (Pr.)

One Size

85696-1


\$130

Other types of vests and belts available, please see Vest and Belts Order Form


Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------

**Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable**

- Turtle Bar System ♦ A: 5", B: 4-8½" Short 98619-25 \$230
  - Turtle Bar System ♦ A: 10", B: 4-8½" Long 98619-26 \$230
- ♦Required with Anatomic Head Support; 5099001, 5099002 or 5099003



- Anatomic Head Support ♦ A:5½ x B:5½ x C:9½" Small 5099001 \$255
  - Head Support Cover (Required) Small 9923-83 \$0
  - Anatomic Head Support ♦ A:6 x B:5½ x C:10" Medium 5099002 \$255
  - Head Support Cover (Required) Medium 9945-83 \$0
  - Anatomic Head Support ♦ A:6½ x B:5½ x C:12½" Large 5099003 \$255
  - Head Support Cover (Required) Large 9967-83 \$0
- ♦Turtle Bar System Short or Long Required; 98619-25 or 98619-26



**XX = Turtle Bar Size**

- 92 = Short Turtle Bar
- 93 = Long Turtle Bar

**Adjustable Headrest without Dynamic Feature - Requires Bracket (see below)**

- Adj. Width Headrest, Large Laterals, Black Cover ♦ 99737xx-85 \$305
  - Adj. Width Headrest, Large Laterals, w/Speakers, Black Cover ♦♦ 99717xx-85 \$405
- ♦♦Custom order; Approx. lead time: 6-8 weeks.  
♦Turtle Bar System Short or Long Required; x98619-25 or x98619-26



9971xxx-xx

9973xxx-xx

**Adjustable Headrest with Dynamic Feature - Requires Bracket (see below)**

- Adj. Width Headrest, Large Laterals, Black Cover ♦ 99736xx-85 \$385
  - Adj. Width Headrest, Large Laterals, w/Speakers, Black Cover ♦♦ 99756xx-85 \$385
- ♦♦Custom order; Approx. lead time: 6-8 weeks.  
♦Turtle Bar System Short or Long Required; x98619-25 or x98619-26



9975xxx-xx

- Headband System for Anatomic Head Support Only One Size 99500 \$230

- Universal Headrest Mount ♦ One Size 87590 \$61

♦Required when using other manufactures Head Support

Allows most Whitmyer, Stealth, Otto Bock headrest, I21 to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation

Other types of headrest available, please see Head Supports Order Form

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------

- Directional Stabilizer Kit ♦♦ One Size 89592-KIT \$215
- ♦♦Custom order; Approx. lead time: 6-8 weeks.

- Tray (Ea.) 1,2 95125 \$300
- Tray (Ea.) 3 95126 \$300



Accessories	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------

- Tray, Grey 1,2 8976350 \$225



- Push Brace, Seat Mounted 1,2 95148-1 \$230
- Push Brace, Seat Mounted 3 95148-3 \$230





Angle Adjustable Push Brace (Ea.) 3 95147 \$210




Helper's Handle, Fold Down (Pr.) One Size 95145 \$80



Helper's Handle, Fixed (Pr.) One Size 95146 \$80


Anti-Tip Device One Size 95953 \$160

Spare Parts	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



Wombat Living Step 1 to 2 Kit, Seat & Back 1 95560 \$280

Wombat Living Step 1 to 2 Kit, Seat & Back 2 95561 \$280



Wombat Living Short to High Back Kit 1 95560-1 \$157

Wombat Living Short to High Back Kit 2 95561-1 \$157

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

**Required, must select a color**

**xx = Color**

Red = 24

Petroleum = 25

Grey = 28

Cover for Standard Back Cushion 1 955611-2xx \$109

Cover for Standard Back Cushion 2 955612-2xx \$109

Cover for Standard Back Cushion 3 955613-2xx \$109

Cover for Contoured Back Cushion, Grey 1 955611-238 \$137

Cover for Contoured Back Cushion, Grey 2 955612-238 \$137

Cover for Contoured Back Cushion, Grey 3 955613-238 \$137

Cover for Seat Cushions 1 955621-2xx \$109

Cover for Seat Cushions 2 955622-2xx \$109

Cover for Seat Cushions 3 955623-2xx \$109

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

Cover w/Foam for Pelvis Support without Sides 1,2 95122-1 \$130

Cover w/Foam for Pelvis Support without Sides 3 95122-2 \$130

Cover w/Foam for Pelvis Support with Sides 1 95123-1 \$154

Cover w/Foam for Pelvis Support with Sides 2 95123-2 \$154

Cover w/Foam for Pelvis Support with Sides 3 95123-3 \$154

Covers for Fixed Side Supports (95121) 1,2 95121-60 \$84

Spare Parts	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



Footplate, Angle & Height Adj. 12" x 7½" 95145-1 \$115

Footplate, Angle & Height Adj. 14" x 9" 85145-2 \$115



Footplate, Angle & Height Adj. ♦ 14" x 9" 95145-2 \$115

Footplate, Angle & Height Adj. ♦ 16" x 9½" 95145-3 \$115

Footplate, Angle & Height Adj. ♦ 18¾" x 9½" 95145-4 \$115

♦ Includes 4" mounting fitting



Footplate with Rod 7½" x 12" 1 95144-1 \$225



<input type="checkbox"/> Footplate with Rod	9" x 14"	2	95144-2	\$225
<input type="checkbox"/> Footplate with Rod	9¼" x 16"	3	95144-3	\$225
<input type="checkbox"/> Footplate with Rod	9½" x 19"	4	95144-4	\$225

---

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

---

<input type="checkbox"/> Batteries for Wombat Living with Power (Set)	One Size	2222012	\$305
◆ <b>Both batteries <u>must</u> be charged at the same time</b>			



<input type="checkbox"/> Charger & Plug for USA/Canada (Ea.)	One Size	2222161	\$320
--	----------	---------	-------

---

Prices in USD & is subject to change without notice

**For Quotes and Orders please send completed Order Form to: [sales.us@etac.com](mailto:sales.us@etac.com) or fax: 1-800-760-0537**