

mustangTM

Exploded views / Part lists

Rev. date (12.2017)

On the following pages you will find the parts used to assemble this product. Please note that some spare parts cannot be purchased individually but are part of a complete kit. Please also note that for some spare parts extra delivery time may occur.



Contents


| | |
|------------------------|---|
| Mustang, standard..... | 3 |
| Accessories | 7 |



| Position | No. | Description | Unit of Measure | Quantity | Ean code |
|----------|----------|--------------------------------|-----------------|----------|---------------|
| 1 | 79438-R | Wheel size 2 mounted, Right | PCS | 2 | 5707292543448 |
| 1 | 79438 | Wheel Ø200, size 2 mounted | PCS | 4 | 5707292527073 |
| 2 | 79436-L | Wheel size 1 Left, Mustang | PCS | 2 | 5707292543417 |
| 3 | 76501-38 | Frame welded (size 1) Yelleow | PCS | 1 | 5707292502834 |
| 3 | 76502-24 | Frame welded (size 2) blue | PCS | 1 | 5707292521545 |
| 3 | 76502-35 | Frame welded (size 2) Antrazit | PCS | 1 | 5707292632852 |
| 3 | 76502-31 | Frame welded (size 2) Red | PCS | 1 | 5707292603319 |
| 3 | 76502-10 | Frame welded (size 2) white | PCS | 1 | 5707292603296 |
| 3 | 76503-27 | Frame welded (size 3), red | PCS | 1 | 5707292548115 |
| 3 | 76503-35 | Frame welded (size 3) Antrazit | PCS | 1 | 5707292632869 |
| 3 | 76503-31 | Frame welded (size 3), Red | PCS | 1 | 5707292603357 |
| 3 | 76503-10 | Frame welded (size 3), White | PCS | 1 | 5707292603326 |
| 3 | 76504-23 | Frame welded (size 4) blue | PCS | 1 | 5707292548122 |
| 3 | 76504-35 | Frame welded (size 4) antrazit | PCS | 1 | 5707292632876 |
| 4 | 76523 | Fix bracket (left) | PCS | 2 | 5707292507570 |
| 5 | 76522 | Fix bracket (right) | PCS | 2 | 5707292507563 |
| 6 | 0701012 | Stop screw M10*12 flat stainl. | PCS | 1 | 5707292533593 |
| 7 | 76715 | Pin Ø8x27mm | PCS | 2 | 5707292533722 |
| 8 | 1145060 | Locking washer, axle stainless | PCS | 4 | 5707292375056 |
| 9 | 2501408 | Black blind ferrule | PCS | 2 | 5707292143921 |
| 9 | 2501410 | Black blind ferrule Ø8,3*12mm | PCS | 2 | 5707292577719 |
| 10 | 0602425 | Insex M6*25, stainless | PCS | 4 | 5707292404596 |
| 11 | 0602205 | Stainless insex M6*16 | PCS | 4 | 5707292404732 |
| 12 | 76720 | Q-release, for frame | PCS | 2 | 5707292533753 |
| 13 | 1902730 | Crosshead screw M5*10, FZB | PCS | 2 | 5707292433084 |
| 14 | 1631060 | Stainless lock nut 6mm | PCS | 6 | 5707292141446 |
| 14 | 1606500 | Lock nut 6 mm | PCS | 6 | 5707292141361 |
| 15 | 1607002 | Lock nut 8mm, A2 | PCS | 3 | 5707292410955 |
| 16 | 1600524 | Stainl. Welding nut w/flange | PCS | 2 | 5707292503541 |
| 17 | 1105080 | Stainless washer 8mm | PCS | 4 | 5707292140364 |
| 20 | 76517 | Frontcolumn (size 1+2) | PCS | 1 | 5707292507525 |
| 21 | 76519 | Bearing shell (size 1+2) | PCS | 2 | 5707292507549 |
| 22 | 76518 | Overlay (size 1+2), frontcolum | PCS | 1 | 5707292507532 |
| 23 | 76531-38 | Stabilisator frame welded siz1 | PCS | 1 | 5707292502957 |
| 23 | 76532-24 | Stabilisator frame welded siz2 | PCS | 1 | 5707292521552 |
| 23 | 76532-35 | Stabilisator frame welded siz2 | PCS | 1 | 5707292632883 |
| 23 | 76532-31 | Stabilisator frame welded siz2 | PCS | 1 | 5707292603623 |
| 23 | 76532-10 | Stabilisator frame welded siz2 | PCS | 1 | 5707292603616 |
| 23 | 76533-27 | Stabilisator frame welded siz3 | PCS | 1 | 5707292548139 |
| 23 | 76533-35 | Stabilisator frame welded siz3 | PCS | 1 | 5707292632890 |
| 23 | 76533-31 | Stabilisator frame welded siz3 | PCS | 1 | 5707292603654 |
| 23 | 76533-10 | Stabilisator frame welded siz3 | PCS | 1 | 5707292603630 |
| 23 | 76534-23 | Stabilisator frame welded siz4 | PCS | 1 | 5707292548146 |
| 23 | 76534-35 | Stabilisator frame welded siz4 | PCS | 1 | 5707292632906 |
| 24 | 76539 | Name plate, Mustang | PCS | 1 | 5707292541833 |
| 25 | 1940007 | Eco syn 4*16 torx | PCS | 8 | 5707292520593 |
| 27 | 6153600 | Knob spring f.Panther/ Mustang | PCS | 1 | 5707292188588 |
| 28 | 0631003 | Insex unders. M10*35, rustfri | PCS | 2 | 5707292539670 |
| 28 | 0614005 | Insex unders.M10*40, rf | PCS | 2 | 5707292424457 |

| | | | | | |
|----|---------|--------------------------------|-----|---|---------------|
| 29 | 76555 | Bush, topcolumn | PCS | 2 | 5707292507846 |
| 29 | 76559 | Bush, size 3+4 topcolumn Ø10 | PCS | 2 | 5707292528537 |
| 30 | 76551 | Top column (size 1+2) | PCS | 1 | 5707292507822 |
| 30 | 76552 | Top column (size 3+4) | PCS | 1 | 5707292518941 |
| 30 | 110227 | Top column (size 3+4 black) | PCS | 1 | 5707292580245 |
| 31 | 76575 | Bush, size 1-2 | PCS | 2 | 5707292520647 |
| 31 | 76553 | Bush, topcolumn size 3-4 | PCS | 2 | 5707292528377 |
| 32 | 76554 | Bush, topcolumn | PCS | 1 | 5707292507839 |
| 32 | 76558 | Bush, topcolumn with M10 | PCS | 1 | 5707292528384 |
| 33 | 107625 | Spindle for angle adj. | PCS | 1 | 5707292571977 |
| 33 | 107634 | Spindle for angle adj. | PCS | 1 | 5707292572004 |
| 33 | 79031 | Spindel w skala for adjustment | PCS | 1 | 5707292401595 |
| 34 | 2200307 | Bloc-o-lift 200n | PCS | 1 | 5707292539656 |
| 34 | 2200308 | Bloc-o-lift 200n | PCS | 1 | 5707292539786 |
| 35 | 3470335 | Release button 4*8mm | PCS | 1 | 5707292146632 |
| 36 | 76706 | Pedal plate size 1+2 | PCS | 1 | 5707292533739 |
| 36 | 115713 | New Pedal plate size 3+4 | PCS | 1 | 5707292606402 |
| 37 | 6118822 | Horizontal release | PCS | 1 | 5707292186744 |
| 38 | 2504500 | Ferrule e20*1,0-2,0 int.black | PCS | 1 | 5707292144126 |
| 40 | 0631002 | Insex unders. M10*30, rustfri | PCS | 2 | 5707292539663 |
| 40 | 0614006 | Insex unders.M10*60, rf | PCS | 2 | 5707292544957 |
| 41 | 0603200 | Insex screw 8*40 | PCS | 1 | 5707292138224 |
| 42 | 0603000 | Insex screw 8*35 | PCS | 2 | 5707292138217 |
| 43 | 1761000 | Counter nut 10*1,0 fine thread | PCS | 1 | 5707292141538 |
| 44 | 1310717 | Nylon washer, nat o17-8,4*2mm | PCS | 1 | 5707292140821 |
| 45 | 1310712 | Nylon washer, nature | PCS | 1 | 5707292356345 |
| 46 | 1310601 | Nylon washer din 125, m8 | PCS | 1 | 5707292557261 |
| 50 | 76700 | Center joint | PCS | 1 | 5707292527127 |
| 51 | 76581 | Centre pole size 1+2 | PCS | 1 | 5707292507983 |
| 51 | 76582 | Centre pole size 3+4 | PCS | 1 | 5707292508010 |
| 52 | 5082922 | Plastbøsning 10/30, Mustang | PCS | 2 | 5707292557766 |
| 52 | 82922 | Bush for Armrest | PCS | 2 | 5707292540539 |
| 53 | 76686 | Flip support with pur size 1-c | PCS | 1 | 5707292508638 |
| 53 | 76687 | Flip support with pur size 2-c | PCS | 1 | 5707292508645 |
| 53 | 76688 | Flip support with pur size 3-c | PCS | 1 | 5707292508652 |
| 54 | 76676 | Side support with pur size 1-a | PCS | 1 | 5707292508522 |
| 54 | 76679 | Side support with pur size 2-b | PCS | 1 | 5707292508560 |
| 54 | 76680 | Side support with pur size 3-b | PCS | 1 | 5707292508577 |
| 55 | 76686 | Flip support with pur size 1-c | PCS | 1 | 5707292508638 |
| 55 | 76689 | Flip support with pur size 2-d | PCS | 1 | 5707292508669 |
| 55 | 76690 | Flip support with pur size 3-d | PCS | 1 | 5707292508676 |
| 56 | 76686 | Flip support with pur size 1-c | PCS | 1 | 5707292508638 |
| 56 | 76689 | Flip support with pur size 2-d | PCS | 1 | 5707292508669 |
| 56 | 76690 | Flip support with pur size 3-d | PCS | 1 | 5707292508676 |
| 57 | 76676 | Side support with pur size 1-a | PCS | 1 | 5707292508522 |
| 57 | 76677 | Side support with pur size 2-a | PCS | 1 | 5707292508539 |
| 57 | 76678 | Side support with pur size 3-a | PCS | 1 | 5707292508553 |
| 58 | 76686 | Flip support with pur size 1-c | PCS | 1 | 5707292508638 |
| 58 | 76687 | Flip support with pur size 2-c | PCS | 1 | 5707292508645 |
| 58 | 76688 | Flip support with pur size 3-c | PCS | 1 | 5707292508652 |

| | | | | | |
|----|---------|--------------------------------|-----|---|---------------|
| 59 | 5540410 | Band Mustang size 1-2-3-4 | PCS | 1 | 5707292527509 |
| 60 | 76596 | Width adjustment size 1 | PCS | 2 | 5707292508171 |
| 60 | 76597 | Width adjustment size 2 | PCS | 2 | 5707292508188 |
| 60 | 76598 | Width adjustment size 3+4 | PCS | 2 | 5707292538918 |
| 61 | 76570 | Extension muff | PCS | 1 | 5707292519023 |
| 62 | 76716 | Front tube, size 1+2 | PCS | 1 | 5707292533760 |
| 62 | 76717 | Front tube, size 3+4 | PCS | 1 | 5707292552310 |
| 63 | 1310550 | Covering stud, black, Gazelle | PCS | 1 | 5707292140722 |
| 64 | 6153700 | Knob spring f.activity centre | PCS | 1 | 5707292188601 |
| 69 | 1130600 | Stainless washer 6mm | PCS | 2 | 5707292140517 |
| 71 | 0701013 | Stop screw M10*16 flat | PCS | 3 | 5707292533609 |
| 72 | 0701012 | Stop screw M10*12 flat stainl. | PCS | 2 | 5707292533593 |
| 73 | 0641621 | Button head w/flange M6*20 FZB | PCS | 2 | 5707292360960 |
| 73 | 0659501 | Button head w/flange M6*16 FZB | PCS | 2 | 5707292356000 |
| 74 | 0602470 | Insex M6*70, stainless | PCS | 4 | 5707292407870 |
| 75 | 1606500 | Lock nut 6 mm | PCS | 4 | 5707292141361 |
| 76 | 0624801 | Insex screw counters. 8*20 | PCS | 4 | 5707292139061 |
| 77 | 1101750 | Washer, conic, M6, FZB | PCS | 4 | 5707292528421 |
| 78 | 1781050 | Top nut M5 / lock rf | PCS | 2 | 5707292530639 |
| 79 | 0658030 | Button head M5*60 stainless | PCS | 2 | 5707292530622 |
| 80 | 76696 | Frontcolumn (size 3+4) | PCS | 1 | 5707292535313 |
| 81 | 76699 | Bearing shell (size 3+4) left | PCS | 1 | 5707292535344 |
| 82 | 76698 | Bearing shell size 3+4 right | PCS | 1 | 5707292535337 |
| 83 | 76697 | Overlay (size 3+4) | PCS | 1 | 5707292535320 |
| 84 | 82950 | Hollow plate 16*32, Wombat | PCS | 1 | 5707292126030 |
| 85 | 1940009 | Eco syn 4*20 torx | PCS | 8 | 5707292550897 |
| 86 | 1940037 | Eco-syn z M4*16 FZB | PCS | 2 | 5707292458650 |
| 86 | 1940017 | Sheet metal screw uhx4,2*13fzb | PCS | 2 | 5707292142597 |
| 87 | 75291 | Hexagon key nv=5 | PCS | 1 | 5707292437655 |
| 88 | 76542 | Bumper, wheel | PCS | 4 | 5707292507709 |
| 89 | 1811205 | Sheet metal screw 4,8*22 phx | PCS | 4 | 5707292530233 |
| 90 | 76718 | Spring f. pedal plate,size 1+2 | PCS | 1 | 5707292583574 |
| 90 | 0310587 | Vridfjeder | PCS | 1 | 5707292556240 |
| 91 | 79436-R | Wheel size 1 Right, Mustang | PCS | 2 | 5707292543424 |
| 92 | 79438-L | Wheel size 2 mounted, Left | PCS | 2 | 5707292543431 |
| 93 | 1310531 | Covering plug,black Ø6,0 mm | PCS | 4 | 5707292396211 |



| Position | No. | Description | Unit of Measure | Quantity | Ean code |
|----------|---------|--------------------------------|-----------------|----------|---------------|
| 100 | 76614 | Saddle with pur size 2 | PCS | 1 | 5707292539298 |
| 100 | 76613 | Saddle with pur | PCS | 1 | 5707292508218 |
| 101 | 3502115 | Tightning handle M6x15mm | PCS | 1 | 5707292146854 |
| 102 | 0380184 | Spec. Locking house M10*1,28mm | PCS | 1 | 5707292366030 |
| 103 | 0701012 | Stop screw M10*12 flat stainl. | PCS | 1 | 5707292533593 |
| 104 | 0506516 | Stainless blind rivet M5*16 | PCS | 4 | 5707292137357 |
| 104 | 0506516 | Stainless blind rivet M5*16 | PCS | 6 | 5707292137357 |
| 105 | 76641 | Seperator plate, leg. Left sz1 | PCS | 1 | 5707292508317 |
| 105 | 76548 | Seperator plate size 2-4 left | PCS | 1 | 5707292518927 |
| 106 | 76642 | Seperator plate, leg sep.right | PCS | 1 | 5707292534552 |
| 106 | 76549 | Seperator plate size 2-4 right | PCS | 1 | 5707292518934 |
| 107 | 76645 | Fix bracket, legseperator | PCS | 1 | 5707292508348 |
| 107 | 118201 | Basic profile, legseperator | PCS | 1 | 5707292613479 |
| 107 | 76561 | Centre joint | PCS | 1 | 5707292507877 |
| 107 | 76564 | Bracket, centre joint | PCS | 1 | 5707292507914 |
| 107 | 76579 | Basic profile, size 4 | PCS | 1 | 5707292543615 |
| 108 | 2507000 | Ferrule 30*20*1,0-2,0 int.blac | PCS | 2 | 5707292144232 |
| 109 | 2504500 | Ferrule e20*1,0-2,0 int.black | PCS | 1 | 5707292144126 |
| 110 | 76545 | Lock tap | PCS | 2 | 5707292507761 |
| 111 | 76544 | Lock cover | PCS | 2 | 5707292507747 |
| 112 | 76536 | Pole (size 1) | PCS | 2 | 5707292507679 |
| 112 | 76608 | Pole (size 2 special) | PCS | 2 | 5707292554970 |
| 112 | 76537 | Pole (size 2) | PCS | 2 | 5707292518910 |
| 112 | 76609 | Pole (size 3) | PCS | 2 | 5707292554956 |
| 112 | 76610 | Pole (size 4) | PCS | 2 | 5707292554963 |
| 113 | 76515 | Torsion spring | PCS | 4 | 5707292507518 |
| 115 | 0310100 | Pressure spring 0,9*8*23*10,5v | PCS | 2 | 5707292135513 |
| 116 | 0602453 | Insex M5*40 stainless | PCS | 2 | 5707292519009 |
| 117 | 0530012 | Cylind.pin Ø8x40 mm. stainless | PCS | 4 | 5707292541024 |
| 119 | 1105060 | Stainless washer 6mm | PCS | 4 | 5707292140357 |
| 120 | 2500516 | Grip 87.00 mm length | PCS | 2 | 5707292356727 |
| 122 | 0603530 | Insex low head M8x16 stainless | PCS | 2 | 5707292412492 |
| 123 | 0609744 | Stainless insex 6*14 ch. | PCS | 4 | 5707292418234 |
| 124 | 0701012 | Stop screw M10*12 flat stainl. | PCS | 2 | 5707292533593 |
| 124 | 0701012 | Stop screw M10*12 flat stainl. | PCS | 1 | 5707292533593 |
| 125 | 76662 | Armrest welded, left | PCS | 1 | 5707292508454 |
| 126 | 76627 | Handle, left | PCS | 1 | 5707292535276 |
| 127 | 76661 | Armrest welded, right | PCS | 1 | 5707292508393 |
| 128 | 76659 | Armrest with pur, right | PCS | 1 | 5707292508386 |
| 129 | 76658 | Armrest with pur, left | PCS | 1 | 5707292508379 |
| 130 | 0701013 | Stop screw M10*16 flat | PCS | 1 | 5707292533609 |
| 131 | 76621 | Armrest holder, armrest | PCS | 1 | 5707292508249 |
| 131 | 76622 | Armrest holder, size 3+4 | PCS | 1 | 5707292552839 |
| 132 | 0602205 | Stainless insex M6*16 | PCS | 12 | 5707292404732 |
| 133 | 76577 | Gatetrainer adapter | PCS | 4 | 5707292535641 |
| 134 | 76578 | Gatetrainer adapter | PCS | 4 | 5707292535665 |
| 135 | 1631060 | Stainless lock nut 6mm | PCS | 12 | 5707292141446 |
| 137 | 72845 | Chest support size 2 | PCS | 1 | 5707292445834 |
| 138 | 1350024 | Distance piece 14*15*8,2mm | PCS | 1 | 5707292140920 |

| | | | | | |
|-----|---------|--------------------------------|-----|---|---------------|
| 139 | 2504000 | Ferrule e15*1,0-2,0 int.black | PCS | 2 | 5707292144102 |
| 139 | 2504000 | Ferrule e15*1,0-2,0 int.black | PCS | 1 | 5707292144102 |
| 140 | 1940016 | Plate screw 4,2*19 FZB | PCS | 8 | 5707292356796 |
| 141 | 73448 | Armholder 1,crocodile | PCS | 2 | 5707292355126 |
| 141 | 73449 | Arm rest plate sz 2,crocodile | PCS | 2 | 5707292355133 |
| 142 | 76626 | Handle, right | PCS | 1 | 5707292508256 |
| 143 | 1607000 | Lock nut 8mm, large pack | PCS | 2 | 5707292141378 |
| 144 | 2503010 | Ferrule d22*1,8-3,5, black | PCS | 2 | 5707292357632 |
| 145 | 3502516 | Tightening handle 45mm M8*16 | PCS | 2 | 5707292146915 |
| 146 | 1750800 | Lock nut, stainless M8 | PCS | 2 | 5707292141484 |
| 147 | 73446 | Retainer forearm support | PCS | 2 | 5707292355102 |
| 148 | 102798 | Armrest adjustment (top) | PCS | 2 | 5707292564559 |
| 149 | 0663041 | Button head M8*40, rf | PCS | 2 | 5707292575388 |
| 150 | 102799 | Armrest adjustment (bottom) | PCS | 2 | 5707292564566 |
| 151 | 0310129 | pressure spring 0,9*11*26mm*7v | PCS | 2 | 5707292438782 |
| 152 | 73467 | Tube for handle Crocodile | PCS | 2 | 5707292356314 |
| 152 | 76588 | Tube part f. armrest | PCS | 2 | 5707292535481 |
| 153 | 2502002 | Ferrule d18*0,8-2,5 black | PCS | 2 | 5707292433947 |
| 154 | 5020910 | Arm rest cover size 1 | SET | 1 | 5707292361158 |
| 154 | 5020920 | Cover for forearm support sz 2 | SET | 1 | 5707292361165 |
| 155 | 128485 | Steel bushing for armrest | PCS | 2 | 5707292795304 |
| 160 | 72830 | Chestsupport size 2 right | PCS | 1 | 5707292421357 |
| 160 | 72840 | Chestsupport size 3 right | PCS | 1 | 5707292421371 |
| 161 | 72831 | Chestsupport size 2 left | PCS | 1 | 5707292421364 |
| 161 | 72841 | Chestsupport size 3 left | PCS | 1 | 5707292421388 |
| 162 | 72738 | Mounting cross, backsupport | PCS | 1 | 5707292422989 |
| 163 | 1102010 | Face plate 8,4*18*2mm, FZB | PCS | 4 | 5707292429766 |
| 164 | 0603528 | Insex lav hoved M8*14 | PCS | 4 | 5707292564474 |
| 165 | 3502516 | Tightening handle 45mm M8*16 | PCS | 1 | 5707292146915 |
| 166 | 2504000 | Ferrule e15*1,0-2,0 int.black | PCS | 3 | 5707292144102 |
| 167 | 116908 | Tube, angle | PCS | 1 | 5707292612595 |
| 170 | 75294 | Mounting cross, back support | PCS | 1 | 5707292509963 |
| 180 | 72850 | Back support | PCS | 1 | 5707292423023 |
| 180 | 5038100 | Upholstered sacrum support 2 | PCS | 1 | 5707292553386 |
| 181 | 76530 | Mounting plate, front support | PCS | 1 | 5707292507617 |
| 182 | 0603534 | Insex low head M8x20 | PCS | 1 | 5707292575302 |
| 183 | 0603532 | Insex lav hoved M8*18 | PCS | 2 | 5707292564481 |
| 184 | 1102010 | Face plate 8,4*18*2mm, FZB | PCS | 2 | 5707292429766 |
| 185 | 76616 | Tube 25x25 to 15x15 adapter | PCS | 1 | 5707292508225 |
| 186 | 116910 | Mounting plate, butt support | PCS | 1 | 5707292612601 |
| 189 | 82894 | Saddle washer | PCS | 1 | 5707292387578 |
| 190 | 0603534 | Insex low head M8x20 | PCS | 1 | 5707292575302 |
| 191 | 6153550 | Knob spring Ø 7.5x5.5 mm | PCS | 1 | 5707292513717 |
| 192 | 0701012 | Stop screw M10*12 flat stainl. | PCS | 1 | 5707292533593 |
| 193 | 1105080 | Stainless washer 8mm | PCS | 1 | 5707292140364 |
| 194 | 0603541 | Insex low head M8*60 | PCS | 1 | 5707292575364 |
| 195 | 76660 | Push brace bracket | PCS | 1 | 5707292534965 |
| 196 | 76646 | Angle joint, guide brace | PCS | 1 | 5707292519528 |
| 197 | 1603500 | Steel nut 8mm, large pack | PCS | 1 | 5707292141286 |
| 198 | 76647 | Brace welded | PCS | 1 | 5707292519535 |

| | | | | | |
|-----|------------|--------------------------------|-----|-----|---------------|
| 199 | 2500520 | Tpe cover Ø20mm | PCS | 1 | 5707292401649 |
| 200 | 2503010 | Ferrule d22*1,8-3,5, black | PCS | 1 | 5707292357632 |
| 210 | 76725 | Armrest holder | PCS | 1 | 5707292534569 |
| 210 | 101009 | Armrest holder, size 2 | PCS | 1 | 5707292552938 |
| 211 | 76727 | Front tube welded, armrest | PCS | 1 | 5707292534583 |
| 211 | 101007 | Front tube welded, size 2 | PCS | 1 | 5707292552976 |
| 212 | 2503010 | Ferrule d22*1,8-3,5, black | PCS | 1 | 5707292357632 |
| 213 | 2500550 | Grip 700mm length for Ø20mm | PCS | 0,5 | 5707292388520 |
| 213 | 2500550 | Grip 700mm length for Ø20mm | PCS | 1 | 5707292388520 |
| 214 | 76510 | Ende prop | PCS | 1 | 5707292502872 |
| 204 | 869560 | Storage pouches | SET | 1 | 5707292538819 |
| 205 | 869810-160 | Cover for chest support, small | SET | 1 | 5707292538789 |
| 205 | 869820-260 | Cover for chest support,medium | SET | 1 | 5707292538796 |
| 210 | 76725 | Armrest holder | PCS | 1 | 5707292534569 |
| 211 | 76727 | Front tube welded, armrest | PCS | 1 | 5707292534583 |
| 212 | 2502500 | Ferrule d20*1,0-2,0 int.black | PCS | 1 | 5707292144027 |
| 213 | 2500550 | Grip 700mm length for Ø20mm | PCS | 0,5 | 5707292388520 |
| 214 | 76510 | End plug | PCS | 1 | 5707292502872 |
| 215 | 869470 | Leg pad | PCS | 1 | 5707292553256 |
| 220 | 869370 | Padded forearm positioner | PCS | 1 | 5707292579676 |